

KeyPrint® 3D Printer Validation

Status: June 24, 2019



We continue to validate printers. Can't find your printer? Call us!
We'll work with you to validate your printer.

Visit our web for future updates!

- Validation complete
- Validation in process

Biocompatible Resins
Available now for sampling

		 KeyModel®	 KeyOrthoModel®	 KeyCast®	 KeyMask™	 KeyGuide™	 KeySplint Soft™ <i>510k-pending</i>
	MSRP	\$199 per kg	\$199 per kg	\$325 per kg	\$175 per 0.5 kg	\$299 per kg	\$425 per kg
 EvoDent E120 (405 nm)	●	●	●	●	●	□	□
 EvoDent S Series (405 nm)	●	●	●	●	●	□	□
 Asiga® MAX (385 nm)	●	●	●	●	●	●	□
 Miicraft 125 (405 nm)	●	●	●	●	●	●	●
 Miicraft Ultra (405 nm)	●	●	●	●	●	●	●
 Carbon® M2 (385 nm)							●
 Structo® DentaForm (405 nm)	□	□	□	□	□	□	□
 Rapid Shape® D30 (385 nm)	□	□	□	□	□	□	□
 Ackuretta Diplo (405 nm)	●	●	●	●	●	●	□
 Bego Varseo® S (405 nm)	●	□	●	□	□	●	□
 Bego Varseo® (405 nm)	●	□	●	□	□	●	□
 MoonRay® S100 (405 nm)	●	□	●	□	□	●	□

*Asiga is a registered trademark of Zydex Pty, Ltd.

*Varseo and Varseo S is a registered trademark of BEGO Bremer Goldschlägerei Wilh. HerbstGmbH & Co. KG

*MoonRay is a registered trademark of SprintRay, Inc.

*Carbon is a registered trademark of Carbon, Inc.

*Structo is a registered trademark of Structo Pte., Ltd.