



KeySplint Hard™

The newest member of the KeySplint family, KeySplint Hard is designed for rigid dental splints. KeySplint Hard provides therapy for bruxism and cases where tooth immobilization is needed.

FEATURES

- BIOCOMPATIBLE MEDICAL DEVICE
- UNITED STATES FDA 510K-CLEARED
- HEALTH CANADA APPROVED
- CE MARKED (CLASS IIA)
- HARD FOR RIGID SPLINTS



BIOCOMPATIBLE



DURABLE



EASY TO POLISH



EASILY CLEANED



ABRASION RESISTANT

#4200004 (0.5kg bottle)

#4220004 (1kg bottle)

*Currently available from multiple dealers. Please contact us

PROPERTY	STANDARD/METHOD	KEYSPLINT HARD
Flexural Strength	ISO 20795-2	60-65 MPa
Flexural Modulus	ISO 20795-2	1510-1600 MPa
Water Sorption	ISO 20795-2	18 µg/mm ³
Elongation at Break	ASTM D638	9%
Biocompatibility	ISO 10993	Pass



keystone
industries

keyprint.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01